

	<b>SenseAMP</b>	<b>SenseAMP Plus</b>	<b>SenseAMP PLUS LMW</b>	<b>RampUP</b>	<b>RampUP Plus</b>
<b>Input RNA</b>	25-250ng Total RNA	25-250ng Total RNA	0.5-30ng LMW RNA	1-25ng Total RNA	1-25ng Total RNA
<b>Rounds of Amplification</b>	1	1	1	2	2
<b>Can amplify intact, FFPE, or prokaryotic RNA</b>	Yes	Yes		Yes	Yes
<b>Can amplify microRNA</b>			Yes		
<b>Strand Orientation of Amplified Product</b>	Sense	Sense	Sense	Sense	Sense
<b>Compatible with Affymetrix 3' expression arrays</b>	Yes, with Genisphere's cDNA Synthesis Kit			Yes	
<b>Compatible with Glass expression arrays</b>	Yes, with Genisphere's Array 900MPX or cDNA Synthesis kits	Yes, with Genisphere's 3DNA kits		Yes, with Genisphere's Array 900MPX or cDNA Synthesis kits	Yes, with Genisphere's 3DNA kits
<b>Compatible with miRNA arrays</b>			Yes, with Genisphere's FlashTag kits		
<b>Compatible with qRT-PCR applications</b>	Yes		Yes	Yes	